

# Dromex



## FIRST AID KITS



FA3



FA3/REFILL



FA7



FA7/REFILL

### Description

Dromex® FA3 and FA7, first aid safety box range features a set of materials and tools that is used for giving emergency treatment to a sick or injured person.

Dromex® FA3 are basic first aid kits, whilst Dromex® FA7 kits are used in incidents where there is blood spill.

These kits are suitable for use indoors, in offices or manufacturing environments.

These kits consist of a portable and wall mountable storage container, that contains the following medical supplies and tools:

Contents of first aid boxes			
Items	Contents	FA3 Refill	FA7 Refill
Item 1	Wound cleaner/antiseptic (100ml)	1	1
Item 2	Swabs for cleaning wounds	1	1
Item 3	Cotton wool for padding (100g)	1	1
Item 4	Sterile gauze (minimum quality 10)	10	10
Item 5	1 pair of forceps (for splinters)	1	1
Item 6	1 pair of scissors (minimum size 100mm)	1	1
Item 7	1 set of safety pins	1	1
Item 8	4 triangular bandages	4	4
Item 9	4 roller bandages (75mm x 5m)	4	4
Item 10	4 roller bandages (100mm x 5m)	4	4
Item 11	1 roll of elastic adhesive (25mm x 3m)	1	1
Item 12	1 Non-allergic adhesive strip (25mm x 3m)	1	1
Item 13	1 Packet of adhesive dressing strips (minimum quantity 10 assorted sizes)	1	1
Item 14	4 First aid dressing (75mm x 100mm)	4	4
Item 15	4 First aid dressings (150mm x 200mm)	4	4
Item 16	2 Straight splints	2	2
Item 17	2 Pairs large & 2 pairs medium disposable latex gloves	2	2
Item 18	2 CPR mouth pieces or similar devices	2	2
Item 19	Biocide 6g sachets		1
Item 20	Absorbent 40g sachets		1
Item 21	Eusol Disinfectant 50ml		1
Item 22	Household Gloves Large		2
Item 23	Household Gloves Medium		2
Item 24	Biohazard Bag		2

Refill kits does not include storage boxes.

### Specifications

Style: White, portable and wall mountable first aid box with contents.

Box material: Steel.

Additional: Replacement contents are available, Dromex code FA3/Refill and FA7/Refill.

Box Dimensions:



### Compliance & Conformity

Conforms to Regulation 7 of the Hazardous Chemical Substances Regulations of the Occupational Health and Safety Act. Act 85 of 1993.

### Packaging, Storage & Obsolescence

- FA-3, sold as 5 units, packaged in a box.
- FA-3/Refill contents sold as 10 units, in a box.
- FA-7, sold as an individual unit, packaged in a box.
- FA-7/Refill contents sold as 10 units, in a box.

### Cleaning & Maintenance

Kit must be checked regularly and restocked if any items are damaged or are out of date.

It is recommended that all kits are kept clean and stored away from water.

### Use

Items individually labelled with instructions where applicable.

### Disposal

All medical waste should be disposed of correctly according to local regulations and good disposal practice. First aid kit contents should be disposed of, in medical waste bins.